



LABFACILITY


TEMPERATURE & PROCESS TECHNOLOGY

Type T Retractable Curly Thermocouple Lead (ANSI)



XE-7034-001

TYPE-T (ANSI)



XE-7034-001

TYPE-T (ANSI)



XE-7034-001

TYPE-T (ANSI)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Available in type T

Colour coded to ANSI

2 Metres in length when fully extended

Terminated with a fitted miniature plug

Polyurethane outer

Specifications

Specifications

Product Code	XE-7034-001
Cable Length	2M fully extended
Cable Type	Polyurethane outer
Connector Type	Thermocouple Plug
Connector Size	Miniature
Colour	Blue
Standards Met	Colour code (cores & jacket) to ANSI MC96.1
Thermocouple Type	T
Termination Type	Miniature Thermocouple Plug