

## Compression and Grub Screw Fitting Types



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

### Bayonet Fittings - Compression and Grub Screw Fitting Types

When fitted, a bayonet fitting spring loads the sensor to keep it in permanent contact with the measured surface. Bayonet fittings are particularly suitable for mounting sensors in plastic extruders, injection moulding machines, etc...

## Specifications

<b>General Description</b>	When fitted, a bayonet fitting spring loads the sensor to keep it in permanent contact with the measured surface. Bayonet fittings are particularly suitable for mounting sensors in plastic extruders, injection moulding machines, etc...
<b>Material</b>	Stainless steel
<b>Cable Entry</b>	Compression & Grub Screw Fitting Types
<b>Probe Diameter</b>	For 3.0mm, 4.5mm and 6.0mm probe diameters