

## Stainless Steel Compression Fittings - Tapered Thread (BSPT)



[Click here to view product](#)

<b>Product Code</b>	XE-5071
<b>General Description</b>	Male in-line stainless steel compression type with nut & olive
<b>Material</b>	Stainless steel, BSPT tapered thread
<b>Thread</b>	1/8"BSPT, 1/4"BSPT, 1/2"BSPT, 3/8"BSPT or 3/4"BSPT available
<b>Probe Diameter</b>	Suitable for both metric and imperial probe diameters from 1.0mm to 12.0mm