

## Brass Compression Fittings - Tapered Thread (BSPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

To suit a wide range of mineral insulated and other probe type sensors offered by Labfacility.

Brass compression fittings with BSPT (British Standard Pipe Taper) threads are widely used for connecting pipes and tubes in various applications.

BSPT threads are tapered, meaning the thread diameter decreases towards the end, enhancing the seal as the fitting is tightened.

Made from brass, these fittings provide good corrosion resistance, durability, and ease of installation, making them suitable for plumbing, HVAC, and low to medium pressure systems involving water, air, and gas.

### Specifications

## Specifications

<b>General Description</b>	Male in-line stainless steel compression type with nut & olive
<b>Material</b>	Brass
<b>Thread</b>	1/8"BSPT, 1/4"BSPT, 3/8"BSPT or 1/2"BSPT available
<b>Probe Diameter</b>	Suitable for both metric and imperial probe diameters from 1.5mm