

## Olives - Stainless Steel



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel olives are a critical component in compression fittings, offering a reliable and durable solution for connecting pipes and tubes in a variety of applications.

Their corrosion resistance and strength make them ideal for demanding environments, ensuring secure, leak-proof connections. Proper installation and regular maintenance are key to their effective performance and longevity.

### Specifications

## Specifications

|                            |  |
|----------------------------|--|
| <b>General Description</b> | Olives   |
| <b>Material</b>            | Stainless-Steel                                    |
| <b>Probe Diameter</b>      | Suitable for probe diameters between 1.0mm & 8.0mm |