



Olives - Brass



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Brass olives are crucial for ensuring secure and leak-proof connections in a variety of applications.

Their corrosion resistance, durability, and ease of installation make them ideal for plumbing, HVAC, industrial, and automotive systems. Proper installation and regular maintenance are key to maintaining their effective performance over time.

Specifications

Specifications

General Description	Olives
Material	Brass
Probe Diameter	Suitable for probe diameters between 1.5mm & 8.0mm