



Locknuts - Brass



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Brass locknuts are vital for ensuring secure, vibration-resistant connections in a variety of applications.

Their corrosion resistance and durability make them suitable for demanding environments.

Selecting the right type and grade of locknut, combined with proper installation and regular maintenance, ensures reliable performance and long-term durability.

Specifications

Specifications

General Description	Locknuts
Material	Brass
Thread	Available in M8x1, 1/8"BSP, 1/4"BSP & 1/2"BSP Thread Types
Sizes	Thickness from 3.0mm to 5.5mm