

## Reducers - Brass



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Brass Reducers are vital components in piping systems, facilitating the transition between different pipe sizes to meet specific application requirements.

Their corrosion resistance, durability, and strength make them suitable for a wide range of environments and applications.

Proper selection, installation, and maintenance of stainless steel reducers ensure efficient and reliable performance.

### Specifications

## Specifications

|                            |  |
|----------------------------|--|
| <b>General Description</b> | A reducer allows for a change in size to meet specific requirements of the application |
| <b>Material</b>            | Brass  |
| <b>Sizes</b>               | BSPP   |