

Solid Drilled Thermopockets



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless Steel Thermopockets of various sizes are designed for use with industrial style temperature probes to provide protection for the sheathing from corrosive media and to facilitate probe replacement without disturbing the process.

From stock, process connections are 1/2 inch BSPP as standard, a 1/8 inch BSPP compression fitting to suit diameter of probe is required to connect sensor to thermopocket.

Specifications

General Description	From stock, process connections are 1/2 inch BSPP as standard, a 1/8 inch BSPP compression fitting to suit diameter of probe is required to connect sensor to thermopocket.
Material	Stainless steel
Type & Stem Diameter	303SS, 12mm Stem
Stem Length	50mm, 100mm, 150mm