


Solid Drilled High Pressure Thermopockets



- NOTES:
1. Solid Bar Stock Drilled Thermowell
 2. Material 316/316L St.Stl
 3. Pocket Rated 5000PSI
 4. Material Certificate Required for each Batch

Dimensions in mm unless specified otherwise

[Click here to view product](#)

Product Code	XE-5439
General Description	Solid Drilled High Pressure Thermopockets
Material	316/316L Stainless steel
Pocket Rated	5000 PSI
Type & Stem Diameter	316SS Tapered (12mm at tip)
Stem Length	100mm, 150mm, 250mm
Standards Met	Material Certificate supplied for each order

