






# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Fine Gauge Bare Wire Exposed Junction Thermocouple (0.025mm conductors) - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ultra-fine Type K welded tip Thermocouple for Optimum Sensitivity.

### Specifications

## Specifications

<b>Product Code</b>	XE-3517-001
<b>General Description</b>	Very fast response to a fraction of a second, Ideal for minimal invasive temperature sensing.
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	200mm
<b>Cable Type</b>	Fine gauge bare wire
<b>Core / Strands</b>	0.001" (0.025mm)
<b>Max. Temperature</b>	315°C