

## Handheld Stainless Steel Air Probe - Type T IEC 4.0mm x 110mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

316 Stainless Steel Air Probe with a vented sheath for general air temperature measurement. Comprising of a handle, 2 metre coiled cable and a mini plug.

### **Need Your Temperature Sensor calibrated?**

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

**Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.**

### **Specifications**

## **Specifications**

<b>Product Code</b>	XE-3801-001
<b>General Description</b>	316 Stainless Steel Air Probe with a vented sheath for general air temperature measurement.
<b>Sensor Type</b>	T IEC
<b>Probe Diameter</b>	4.0mm
<b>Probe Length</b>	110mm
<b>Sheath Material</b>	316 Stainless Steel
<b>Cable Length</b>	2m Coiled Cable
<b>Termination Type</b>	Fitted mini plug
<b>Max. Temperature</b>	+400°C