

Handheld Stainless Steel Air Probe - Type T IEC 4.0mm x 110mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

316 Stainless Steel Air Probe with a vented sheath for general air temperature measurement. Comprising of a handle, 2 metre coiled cable and a mini plug.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications

Specifications

Product Code	XE-3801-001
General Description	316 Stainless Steel Air Probe with a vented sheath for general air temperature measurement.
Sensor Type	T IEC
Probe Diameter	4.0mm
Probe Length	110mm
Sheath Material	316 Stainless Steel
Cable Length	2m Coiled Cable
Termination Type	Fitted mini plug
Max. Temperature	+400°C