

## **Handheld Stainless Steel Penetration Probe - Type K JIS 1.6mm x 60mm**

Image  
Coming  
Soon

316 Stainless Steel Penetration Probe with a pointed tip for liquid and semi-solid temperature measurement .  
Comprising of a handle, 2 metre coiled cable and a mini plug.

### **Need Your Temperature Sensor calibrated?**

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

**Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.**

### **Specifications**

## **Specifications**

<b>Product Code</b>	XE-3879-001
<b>General Description</b>	Penetration Probe with a pointed tip
<b>Sensor Type</b>	K JIS
<b>Probe Diameter</b>	1.6mm
<b>Probe Length</b>	60mm
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+400°C