

Handheld Stainless Steel Penetration Probe - Type T JIS 3.3mm x 300mm

Image
Coming
Soon

316 Stainless Steel Penetration Probe with a pointed tip for liquid and semi-solid temperature measurement .
Comprising of a handle, 2 metre coiled cable and a mini plug.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications

Specifications

Product Code	XE-3882-001
General Description	Penetration Probe with a pointed tip
Thermocouple Type	T JIS
Probe Diameter	3.3mm
Probe Length	300mm
Cable Length	2m
Cable Termination	Miniature Plug
Max. Temperature	+400°C