

## Surface Temperature Probe - Type K IEC - 90mm long with 10mm tip



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug.

## Specifications

<b>Product Code</b>	XE-3818-001
<b>General Description</b>	Surface temperature fast response thermocouple with ceramic tip and coiled element
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	+900°C
<b>Dimensions</b>	Probe O6 x 90mm long stem with O12 x 10mm tip. Coiled element diameter 7mm