



# IEC Fine Wire Disc Thermocouple, PTFE - Type K

XE-3515-001  
**FINE WIRE DISC THERMOCOUPLE  
TYPE K**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Copper disc thermocouple with high temperature adhesive kapton strips (13x18mm x 5) for quick and easy attachment to surface/process.

## Specifications

<b>General Description</b>	Fast response – ideal for surface measurement
<b>Sensor Type</b>	Copper disc, 6.35mm diameter x 0.3mm thick
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	1m
<b>Cable Type</b>	PTFE insulated twisted wire
<b>Core / Strands</b>	1/0.2mm
<b>Termination Type</b>	Bare tails
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-0°C
<b>Standards Met</b>	IEC-584-3