


## Glassfibre insulated ANSI Exposed Junction Thermocouple with Miniature Plug - Types K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Need Your Temperature Sensor calibrated?**

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

**Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.**

**Specifications****Specifications**

<b>General Description</b>	Fast response, welded exposed junction
<b>Sensor Type</b>	Good temperature resistance (not suitable in fluids)
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1, 2, 5 or 10 metre
<b>Cable Type</b>	glassfibre insulated flat pair cable
<b>Core / Strands</b>	1/0.3mm conductors
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>Tolerance</b>	ANSI MC96.1 Class 2
<b>Standards Met</b>	ANSI MC96.1