

Mineral Insulated Thermocouple with Miniature or Standard Plug ANSI - Types K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

Thermocouple Type	K
Thermocouple Junction	Hot junction insulated from sheath
Probe Diameter	1.0, 1.5 & 3.0 mm
Probe Length	Various
Sheath Material	Stainless steel 310 sheath
Termination Type	Miniature or Standard Plug
Max. Temperature	+1100°C (+750°C for 1.0mm)
Standards Met	ANSI