

Mineral Insulated Thermocouple with Miniature or Standard Plug ANSI - Types K



[Click here to view product](#)

Product Code	XE-3010
Thermocouple Type	K
Thermocouple Junction	Hot junction insulated from sheath
Probe Diameter	1.0, 1.5 & 3.0 mm
Probe Length	Various
Sheath Material	Stainless steel 310 sheath
Termination Type	Miniature or Standard Plug
Max. Temperature	+1100°C (+750°C for 1.0mm)
Standards Met	ANSI