

## Pt100 4 wire class B Resistance Thermometer, KNE Head



[Click here to view product](#)

<b>Product Code</b>	XE-3652
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
<b>Probe Diameter</b>	6mm
<b>Probe Length</b>	200mm, 400mm or 600mm
<b>Sheath Material</b>	316 stainless steel
<b>Head Type</b>	KNE
<b>Termination Type</b>	4 wire connection block, M20 x 1.5mm cable entry (gland included)
<b>Max. Temperature</b>	+450°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP67

<b>Accuracy</b>	Class B
<b>Number of Wires</b>	4

