

Handheld Mineral Insulated General Purpose Probe - Type K,T,PRT



[Click here to view product](#)

| | |
|----------------------------|--|
| Product Code | XE-3806 |
| General Description | Rounded tip suitable for semi solid and liquid temperature measurement |
| Thermocouple Type | K or T |
| Calibration | IEC, ANSI or JIS Calibration or PRT (Black) |
| Probe Diameter | 1.5mm, 3.0mm & 6.0mm |
| Probe Length | 100mm, 120mm & 300mm |
| Cable Length | 2m Coiled Cable |
| Termination Type | Fitted Mini Plug |
| Max. Temperature | +850°C |