


Handheld Stainless Steel Penetration Probes - Type K,T,PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

316 Stainless Steel Penetration Probe with a pointed tip for liquid and semi-solid temperature measurement .
Comprising of a handle, 2 metre coiled cable and a mini plug.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications

General Description	Penetration Probe with a pointed tip
Sensor Type	K, T or PRT
Calibration	IEC, ANSI or JIS
Probe Diameter	1.6mm, 3.3mm or 6.0mm
Probe Length	60mm, 120mm or 300mm
Cable Length	2m
Cable Termination	Miniature Plug
Max. Temperature	+400°C