


100 Ohm Wire Wound Detectors



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Please note: Class 1/10 DIN lead times will vary - please enquire or you will be contacted with a lead time after placing your order.

Specifications

Specifications

General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	+650°C
Min. Temperature	-200°C
Tolerance	Class A, Class B or 1/10 DIN. 1/10DIN relative to IEC 60751 over the range 0°C to 150°C (1/10 of Class B)
Detector Dimensions	Various width / length combination available