






LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Pipe Surface Jubilee Clip Style Thermocouple, Glassfibre stainless steel overbraid - Type K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Stainless steel pipe clamp sensor, jubilee clip style, with adjustable diameter
Thermocouple Type	K or J
Thermocouple Junction	Grounded tip for fast response
Cable Length	2m
Cable Type	Glassfibre stainless steel braided extension cable
Core / Strands	7/0.2mm
Termination Type	Bare tails
Max. Temperature	+350°C
To Fit Pipe Diameters	13-20mm, 22-30mm, 30-45mm and 40-55mm