

IEC Hypodermic Tip Thermocouple, PTFE - Type T



[Click here to view product](#)

Product Code	XE-9511-001
Sensor Type	chisel point Thermocouple
Thermocouple Type	T
Thermocouple Junction	Grounded hot junction
Probe Diameter	1.6mm
Probe Length	100mm
Sheath Material	316 stainless steel sheath
Cable Length	1.5m
Cable Type	PTFE insulated cable
Core / Strands	7/0.2mm

Termination Type	Bare wire tails (+Brown/-White)
Max. Temperature	+250°C
Min. Temperature	-75°C