


IEC Autoclave Thermocouple - Duplex Type T



[Click here to view product](#)

Product Code	XE-9510-001
General Description	First 2.5m section in a flexible stainless steel conduit
Thermocouple Type	Type T Duplex
Thermocouple Junction	2 x insulated tip junctions
Probe Diameter	4mm
Probe Length	250mm
Sheath Material	Stainless steel (316) sheath
Cable Type	Cable is silicon rubber throughout
Max. Temperature	+200Å°C
Min. Temperature	-50Å°C