

## Probe Disinfectant Wipes (Brand may vary). Pack of 200



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Please note: Packaging / Brand may vary - specifications will be the same.**

Probe Disinfectant Wipes are effective against 99.9% of gram-positive and gram-negative bacteria, as well as most viruses, fungi and other micro-organisms including spores. They are pre-dosed with a dual active formulation ideal for providing effective infection control and preventing cross contamination of probes and utensils.

**Areas of use:**

Designed for the disinfection of small surfaces, Pal Probe Disinfectant wipes are ideal for use in food environments on temperature probes and small utensils. They are also suitable for use in laboratories and clinical environments.

Pocket sized tub (Probe sizes: 130 x 125mm)

Quantity per pack - 200

Tub size - 85mm dia x 170mm height.

Compliant

Controlled dosage ensures effective disinfection for every use.

Reduces the risk of cross contamination.

Blue colour wipe is easily visible in food preparation environments.

Convenient

Replaces the need to continuously wash probes and utensils.

Handy pocket pack fits in chef's tunic or utensil drawer.

Fast to use and easy to use.

Cost effective

Saves valuable time and reduces the need for multiple probes and utensils.

(Use biocides safely. Always read the label and product information sheet before use.)

## Specifications

<b>Product Code</b>	XE-6500-001
---------------------	-------------