



**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## Certificate of Conformity

# ISO9001

**CERTIFICATE OF  
CONFORMITY**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A document certified by a competent authority that the supplied good or service meets the required specifications.

Also called certificate of conformance or certificate of compliance.

## Specifications

<b>Product Code</b>	PO-COC
---------------------	--------