

## Mineral Insulated Pt100 Probe With Extension Lead



[Click here to view product](#)

<b>Product Code</b>	XE-3480
<b>General Description</b>	High temperature, rugged, flexible Pt100 probe
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
<b>Probe Diameter</b>	3 or 6mm
<b>Probe Length</b>	150 or 250mm
<b>Sheath Material</b>	Mineral insulated probe, 316 stainless steel sheath
<b>Pot Type</b>	Plain
<b>Cable Length</b>	1m
<b>Cable Type</b>	Teflon® insulated, screened lead
<b>Core / Strands</b>	7/0.2mm

<b>Cable Termination</b>	4-wire tails, allows connection to any PT100 instrument 2, 3 or 4 wire
<b>Max. Temperature</b>	+500°C
<b>Min. Temperature</b>	-50°C
<b>Accuracy</b>	Class B