



## Ceramic Sheathed Pt100 Inserts



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

|                            |  |
|----------------------------|--|
| <b>General Description</b> | Pt100 thin film element housed in ceramic sheath |
| <b>Element</b>             | Pt100  |
| <b>Cable Length</b>        | 50mm or 450mm                                    |
| <b>Cable Type</b>          | PTFE insulated wires                             |
| <b>Core / Strands</b>      | 7/0.15mm   |
| <b>Cable Termination</b>   | Bare Tails                                       |
| <b>Max. Temperature</b>    | +250°C   |
| <b>Min. Temperature</b>    | -50°C  |
| <b>Tolerance</b>           | Class B  |
| <b>Detector Dimensions</b> | Ceramic 5mm dia x 35mm long                      |
| <b>Number of Wires</b>     | 2 or 4   |