



## Ceramic Sheathed Pt100 Inserts - 4-Wire with 450mm Lead Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

**Specifications**

<b>Product Code</b>	XE-3667-001
<b>General Description</b>	Pt100 thin film element housed in ceramic sheath
<b>Element</b>	Pt100
<b>Cable Length</b>	450mm
<b>Cable Type</b>	PTFE insulated wires
<b>Core / Strands</b>	7/0.15mm
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class B
<b>Detector Dimensions</b>	Ceramic 5mm dia x 35mm long
<b>Number of Wires</b>	4