



# Hygienic Pt100 Resistance Thermometer, 1.5" RJT Style Fitting without Transmitter - 4-wire, 6.0mm x 75mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

RJT unions (Ring Joint Type) are often referred to as a standard union. The compromise of an o-ring style gasket which leaves a small crevice where the liner and male part of the union overlap.

Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.

### Specifications

## Specifications

<b>Product Code</b>	XE-5523-001
<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
<b>Sensor Type</b>	Pt100, 4 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 4 wire single element class A
<b>Transmitter</b>	None
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm