


Hygienic Pt100 Resistance Thermometer, 1.5" RJT Style Fitting with fitted Transmitter - 3-wire, 6.0mm x 125mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

RJT unions (Ring Joint Type) are often referred to as a standard union. The compromise of an o-ring style gasket which leaves a small crevice where the liner and male part of the union overlap.

Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.

Specifications

Specifications

Product Code	XE-5526-001
General Description	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
Sensor Type	Pt100, 3 wire single element class A
Head Type	KNEF
Connection Type	Pt100, 3 wire single element class A
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Support Dia (D2)	8.0mm
Support Length (L2)	50mm