



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Hygienic Pt100 Resistance Thermometer, 1.5" Tri-Clamp fitting with Fitted Transmitter - 3-Wire, 6.0mm x 75mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.

Specifications

Specifications

Product Code	XE-5529-001
Sensor Type	Pt100, 3 wire single element class A
Probe Diameter	6.0mm
Probe Length	75mm
Head Type	KNEF
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Number of Wires	3