




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Hygienic Pt100 Resistance Thermometer, 1.5" Tri-Clamp fitting with Fitted Transmitter - 3-Wire, 6.0mm x 75mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.

### Specifications

## Specifications

<b>Product Code</b>	XE-5529-001
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	75mm
<b>Head Type</b>	KNEF
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Number of Wires</b>	3