


Pt100 3 wire class B Resistance Thermometer, KNE Head with Transmitter



[Click here to view product](#)

| | |
|----------------------------|--|
| Product Code | XE-5506 |
| General Description | Pt100 probe with KNE terminal head fitted with a transmitter |
| Sensor Type | Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 3 wire |
| Process Connection | 1/2" BSP x 6mm |
| Calibration | IEC |
| Probe Diameter | 6mm |
| Probe Length | 150mm |
| Sheath Material | 316 stainless steel |
| Head Type | KNE |
| Max. Temperature | +450°C (connection head @ 170°C) |

| | |
|-------------------------|---|
| Min. Temperature | -50°C |
| IP Rating | IP67 |
| Accuracy | Class B |
| Transmitter | Yes |
| Range | Various ranges of transmitter available |
| Number of Wires | 3 |