

Pt100 4 wire class B Resistance Thermometer, Compact KNS Head



[Click here to view product](#)

Product Code	XE-5510
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Process Connection	1/2" BSP x 6mm
Probe Diameter	6mm
Probe Length	Various probe lengths
Sheath Material	316 stainless steel
Head Type	KNS
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67

Accuracy	Class B
Number of Wires	4