





LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PFA Insulated ANSI Fine Gauge Exposed Junction Thermocouple (0.076mm conductors) - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

Product Code	XE-2006-001
General Description	Fine Gauge Exposed Welded Tip Thermocouple
Sensor Type	Fast response, welded tip exposed junction
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	1m
Cable Type	Teflon® PFA insulated wires
Core / Strands	0.003" (0.076mm) conductors
Cable Termination	Bare tails
Max. Temperature	+260°C
Min. Temperature	-75°C
Standards Met	ANSI colour coded