

## Brass Compression Fittings - NPT Thread (NPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

To suit a wide range of mineral insulated and other probe type sensors offered by Labfacility.

Brass compression fittings with NPT (National Pipe Thread) threads are commonly used for connecting pipes and tubes in various applications.

NPT threads are tapered, which creates a tighter seal as the fitting is tightened.

Brass, known for its corrosion resistance, durability, and machinability, makes these fittings suitable for plumbing, HVAC, and low to medium pressure applications involving water, air, and gas.

### **Specifications**

## **Specifications**

<b>General Description</b>	Male in-line compression type with nut & olive
<b>Material</b>	Brass
<b>Thread</b>	1/8" or 1/4" NPT available
<b>Probe Diameter</b>	To suite 3.0mm or 6.0mm probe diameters