



## 4 Core Silicone Rubber Insulated, Rounded Section, PRT Sensor Cable / Wire



25M



**4 CORE**  
XE-2798-025

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**4 Core Silicone Rubber Insulated, Rounded Section, PRT Sensor Cable / Wire**

Can be used to extend up to 4 wire RTD configurations. Commonly used with Pt100? thermometers but suitable with other types Pt130/500/1000? etc.

**Specifications****Specifications**

<b>General Description</b>	4 core silicone rubber insulated 7/0.2mm tin plated copper conductors, twisted together with silicone rubber outer jacket
<b>Cable Length</b>	10, 25m
<b>Cable Type</b>	Silicone rubber insulated
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	Silicone rubber insulated, tin plated copper conductors, twisted together with silicone rubber outer jacket
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-40°C
<b>Number of Wires</b>	4
<b>Standards Met</b>	Cores colour coded 2 x Red / 2 x White to IEC-751