



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## GM1311 High Performance Professional Digital Food Thermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PLEASE NOTE: Lead time will be 7-10 working days when choosing calibration**

GM1311 High Performance Professional Digital Food/BBQ/Meat Thermometer

Environment temperature :0 to 40°C (32 to 104°F)

Environment humidity:10-90%RH

Rate Voltage: 3V

1 x Food Thermometer

1 x English user manual

1 x CR2032

### Specifications

## Specifications

<b>Product Code</b>	XE-4002-001
<b>Accuracy</b>	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
<b>Resolution</b>	0.1°C / °F
<b>Colour</b>	White
<b>Dimensions</b>	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
<b>Function #1</b>	LCD display
<b>Function #2</b>	Temperature lock
<b>Function #3</b>	Max/Min measuring
<b>Function #4</b>	Conversion between Celsius degree and Fahrenheit degree
<b>Function #5</b>	Auto power off
<b>Max. Temperature</b>	300°C (572°F)
<b>Min. Temperature</b>	-50°C (-58°F)
<b>Probe Length</b>	149.5 x 6mm



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Weight</b>	Weight: 0.95 oz (27 g)
<b>Power Supply</b>	CR2032 button battery x 1 (included)