



## IEC Type K Magnetic Strip Thermocouple, PTFE Insulated with Miniature Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A flexible magnetic strip thermocouple capable of being attached to most ferro-magnetic surfaces, the large area of the magnet ensures good adhesion allowing simplicity of installation. Useful in a wide range of temperature monitoring applications, particularly in test and development.

## Specifications

<b>General Description</b>	Ideal for ferrous metal surface measurement
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Magnet</b>	50 x 25 x 1mm magnetic strip
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	PTFE insulated wires, twin twisted
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584
<b>Max. Temperature</b>	100°C
<b>Min. Temperature</b>	-50°C