




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PFA Insulated ANSI Exposed Junction thermocouple with Moulded-on Plug - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

General Description	Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids, gasses, steam etc.
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	1, 2 or 3m
Cable Type	PFA Teflon® twin twisted pair
Core / Strands	1/0.3mm
Termination Type	moulded-on type K ANSI miniature thermocouple plug
Max. Temperature	260°C
Min. Temperature	-75°C
Tolerance	Tolerance Class 1 for optimum accuracy
Standards Met	ANSI MC96.1 Colour Coded