




IEC Adjustable Bayonet Thermocouple, Glassfibre Stainless Steel Overbraided Cable - Type K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Adjustable Bayonet Thermocouple Type K or J IEC & Type J BS

Adjustable Bayonet Thermocouple with 200mm spring with adjustable fixing cap (11.5mm I.D. single slot).

Specifications
Specifications

General Description	Thermocouple with adjustable insertion depth
Sensor Type	6.0mm diameter tip
Thermocouple Type	K, J IEC J BS
Thermocouple Junction	Grounded tip for fast response
Cable Length	2, 3m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Cable Termination	Bare wire or standard plug
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	IEC-584-3 / BS
Bayonet cap	11.5mm I.D. single slot
Connector Type	Plug
Connector Size	Standard
Thread	M6 x 1mm thread
Spring	200mm