

## Master Class Folding Cooking Thermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PLEASE NOTE: Lead time will be 7-10 working days when choosing calibration**

Perfect for precise cooking of beef, veal, lamb, pork, chicken and turkey, ensure you have a delicious roast joint of meat for your Sunday lunch, impressing both friends and family.

The easy to read LCD display clearly shows the internal temperature measured by the stainless steel probe, with a temperature range from -45°C to 200°C and -49°F to 392°F.

This folding thermometer is also great for easily storing in kitchen draws, with the folding mechanism reducing the overall size.

### Specifications

## Specifications

<b>Product Code</b>	XE-3966-001
<b>General Description</b>	The easy to read LCD display clearly shows the internal temperature measured by the stainless steel probe
<b>Sensor Type</b>	Folding Stainless Steel Probe
<b>Max. Temperature</b>	200°C (392°F)
<b>Min. Temperature</b>	-45°C (-49°F)
<b>Dimensions</b>	25.8 x 8.8 x 2.2 cm
<b>Power Supply</b>	1 x LR44 battery (included)
<b>Weight</b>	100g