


## Air Probe with IEC Moulded-on Miniature Plug - Types K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Air Probe with Moulded-on Miniature type K or T IEC Plug**

Moulded plug gives a robust construction as well as providing a waterproof termination.

**Specifications****Specifications**

<b>General Description</b>	vented sheath for general air temperature measurement
<b>Sensor Type</b>	Mineral insulated
<b>Thermocouple Type</b>	K, T
<b>Probe Diameter</b>	4.0mm
<b>Probe Length</b>	110mm
<b>Termination Type</b>	Miniature moulded flat pin plug, colour coded in accordance IEC specification
<b>Max. Temperature</b>	+250°C (Moulded-on plug +110°C)
<b>Standards Met</b>	IEC
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature (moulded-on)
<b>Probe Material</b>	vented sheath, 316 stainless steel