



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Ceramic Sheathed Pt1000 Grade A Inserts - 2-Wire with 50mm Cable Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3668-001
General Description	Pt1000 Grade A thin film element housed in ceramic sheath
Element	Pt1000
Cable Length	50mm
Cable Type	PTFE insulated wires
Core / Strands	7/0.2mm
Max. Temperature	+250°C
Min. Temperature	-50°C
Detector Dimensions	Ceramic 5mm dia x 50mm long
Number of Wires	2