




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Ceramic Sheathed Pt1000 Grade A Inserts - 2-Wire with 50mm Cable Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

**Specifications**

<b>Product Code</b>	XE-3668-001
<b>General Description</b>	Pt1000 Grade A thin film element housed in ceramic sheath
<b>Element</b>	Pt1000
<b>Cable Length</b>	50mm
<b>Cable Type</b>	PTFE insulated wires
<b>Core / Strands</b>	7/0.2mm
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Detector Dimensions</b>	Ceramic 5mm dia x 50mm long
<b>Number of Wires</b>	2