




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Type K 9kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature or Standard Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Rapid, reliable temperature sensing on ferrous surfaces
Sensor Type	Type K thermocouple to IEC 584
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Termination	1.5 or 3 metres PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
Termination Type	Fitted Mini or Standard Plug
Max. Temperature	+ 250°C
Min. Temperature	-50°C